


Activity

EpiPen® Directions

Skills: Text Use, Numeracy, Document Use
Profiled: Counsellor
Other Applications: All

Important Letter from a Parent - Attachment


How to Use the EpiPen® Auto-Injector



The **EpiPen® Auto Injector** is a spring-activated drug delivery system with a hidden needle. This disposable injector delivers one dose of medication. A person should carry more than one EpiPen® Auto-Injector in case more than one dose is needed.

To Use the EpiPen® Auto-Injector

1. Hold unit with black tip pointing downward
2. Remove grey safety cap
3. Place black tip against thigh at 90°
4. Press hard into thigh until auto-injector clicks
5. Hold for 10 seconds
6. Massage area for several seconds
7. Seek medical attention immediately (911)



After use - bend needle against something. Never touch it.
Put in plastic case and take to hospital with patient.